











A NEW ERA OF **DIGITAL TELECARE**

ALARM IS SENT WHEN

Entry room Fall Detected Leaving bed

Too long in bathroom Loud noise (shout for help) Leaving chair

Planned night visits Getting into bed Alarms when offline/online



Furniture scanning gives freedom to use rooms as you wish.

Anonymization the image is annonomised to protect the dignity of the individual Two-way speech, It is possible to have a good conversation between the user and the monitoring centre via the audio unit. The audio unit also notifies the room when the camera is paused.

Simple installation, Mounted on the wall with the included wall bracket, connect to the power supply and it is ready for use. Built-in 4G and support for WiFi and wired networks Disable and enable Nightowl. The night owl can be shut down using a QR code that indicates how long the owl will be inactive.

Night vision Integrated IR light provides an equally good camera image at night, as in daylight.

YOUR SAFETY

AT HOME AND FALLS SOLUTION

STAY AT HOME WITH PEACE OF MIND

- Hard falls
- Slow falls / slide & slipping
- •Fall out of a chair





- Out of Bed
- Sit up in Bed





- •Leave chair
- Exit room
- Too long in bathroom





Loud noise / Call for Help





- Planned digital
- Supervision visits






Alarms when offline/online





AI PHERIPHERALS

When we upgrade POE power supply to 25W, we will re-arrange USB-C to be the preffered audio connection. This will make the USB 2.0 available for a standard Climax RF dongle supporting Pull Cord, **Door Sesnors and Wrist Alarms** that could bring added value to Nattugla.



Climax RF **(0,45W)** Dongle enable any devices with USB serial port to be intergrated into the network with the dongle via virtual COM port, USB to UART interface, providing easy compatibility solution for existing products.

4G DIGITAL LIFELINE WITH THE PENDANT

FEATURES

Mobility Check, Automatic Check-in report, "Help Arrived" Confirmation, Fault Reporting. Silent emergency call.

BATTERY

Low battery monitoring. Battery backup in case of power failure 48 Hrs.

TESTING

Test button for easy access to trigger test report.



COMMUNICATION

Two way hands-free voice communication. Long communication range.

Integrate multiple devices/sensors: Wrist Transmitter, Emergency Pendant, Fall Sensor, PIR Motion Sensor, Panic Button, Smoke Detector, Water Leakage Sensor, and CO Sensors, etc.

YOUR SAFETY **OUTSIDE**

ENJOY EVERYDAY LIFE, IN COMPLETE SAFETY

1 ONE PRESS TO RAISE THE ALARM

In case of need and 24 hours a day, a simple press on the red button of the watch is enough to launch the alarm. This is effective in all situations, such as in the event of a fall outside, discomfort, aggression or feelings of anxiety. The watch is completely autonomous and does not require a case or a mobile phone.

2 A FEW SECONDS TO CONNECT

We talk to you via the watch's built-in speaker and microphone. If necessary, we can geo-locate your GPS position to provide emergency services with accurate information. If the situation requires it. We call a relative, a neighbour or the emergency services to come to your aid.



THE ANYWHERE WATCH



ALARM FUNCTION

One press of the red button to send the alarm



36H OF BATTERY LIFE

We advise you to charge it every night



LIGHT TO WEAR

The watch is discreet and easy to wear



GEOLOCATION

Relatives can receive the GPS position by SMS



இ EASY CHARGING

A special senior magnetic charger is provided



TIME FUNCTION

And of course, like any watch, it shows the time

THE ANYWHERE GPS TRACKER Features:

SOS Alarms

If the SOS button is pressed an alarm is created to our 24/7 monitoring centre Battery 36/48 Hours

Real-time location monitoring

The device can use its data connection to automatically update its alarm status and location in the online tracking platform.



The pebbell has a waterproof classification of IP-66 which means it even works when dropped in a glass of water.



The device can receive incoming calls and initiate emergency calls by using the call button





HOW TO SUBSCRIBE?



Call us to subscribe by phone.

Our advisors will answer you from Monday

